

PURION Service KIT

...is characterized by convienient handling and optimal support of the disinfection performance. The durability of the PURION products is further improved.



manufacturer	PURION® GmbH
type	PURION Service KIT
Cleaning tissue	4 x 100% cotton
Cleaning solution	200 ml. citric acid (50%)
Safety gloves	4 pair - nitril rubber gloves chemical resistance

The PURION Service KIT enables a reliable and convenient cleaning of the immersion tube system. Cleaning should be carried out latest when the UV lamp is replaced at the end of the lifetime or if operation conditions requires it. ¹

The special cleaning tissues made from 100% cotton guarantee that the sensitive quartz glass is not damaged by scratches during the cleaning procedure.

The cleaning solution ensures that after the cleaning no residues of soap can lead to foam formation. Therefore the PURION Service KIT is suitable to drinking water applications.

The chemical resistant Nitril-gloves provide optimal protection during cleaning.

Disposal can be carried out as household waste.4

Application of the PURION Service KIT:3

Product line PURION 400-2501
Product line AIRPURION

Product line UV Set PURION

Advantages

- no damage of the immersion tube system via scratches
- suitable for drinking water applications (no foam formation due to soap residues)
- cost efficiency

info@purion.eu · www.purion.eu

¹ See instruction manual of the UV-plant

² This risk is present with common household cleaners.

³ For all configurations.

⁴ Before disposal of the cleaning solution sufficient thinning is to be ensured (1:1)